



**Kerala State Drugs and Pharmaceuticals Ltd**  
**(A Govt. of Kerala Enterprise)**  
**KALAVOOR, ALAPPUZHA**

Phone :0477-2258184 (Extn:207) ; Email :purchase@ksdp.in  
Web :http://www.ksdp.co.in

**KSDP/PS/EQ/2025-26/400274/T-138**

**17-12-2025**

**Notice Re Inviting Email Quotation for the Supply of Sticker label**  
**Gluteraldehyde Solution 5 ltr Can (In sheet) IP 2% KMSCL – 6600 Nos**

Quotations are invited for the supply of undermentioned goods/ services as per the attached specifications on FOR destination basis at our factory site at Kalavoor,Alapuzha, Kerala state

SI No.	ITEM CODE & DESCRIPTION	UNIT	QUANTITY	EMD	TENDER FEES
1	41400071 Notice Inviting Email Quotation for the Supply of Sticker label Gluteraldehyde Solution 5 ltr Can (In sheet) IP 2% KMSCL – 6600 Nos	NOS	6,600 NOS	♦ Not Applicable	♦ Not Applicable

Quotations should be submitted as per the Proforma given below on your letter head.

SI No.	Name of item	Make	Rate per unit (including freight, if any)	GST %	Offer validity	Remarks, if any
--------	--------------	------	---	-------	----------------	-----------------

**Due Date/Time : 20-Dec-25 / 1 PM**

**Opening Date/Time: 20-Dec-25 / 1.30 PM**

**Offer validity\*:** Minimum 7 days offer validity from the date of closure of bid submission.

Quotation received with offer validity less than 7 days from the date of closure of bid submission will be entirely rejected.

**Please send your lowest offers of the item to our e-mail ksdpender@gmail.com before 1 PM, 20/12 /2025.**

The Quotation should be submitted through a password protected excel sheet. Please share your Password to our email ksdpender@gmail.com@1.30 PM, 20.12.2025.

**TERMS & CONDITIONS**

- 1.Payment terms :-** 30 days after the receipt of the material along with documents, subject to QC approval.
- 2.Mode of payment:-** E-Payment.
- 3.Delivery Period:-** Supply should be effected within 15 days on award of PO
- 4.Special instructions:-** Quotation number Should be mentioned in the subject line.

**HOD - Purchase**

170 mm

100 mm

5 Litre

## GLUTARALDEHYDE SOLUTION I.P. 2%

Composition: Glutaraldehyde 2%w/v  
Purified Water I.P. : q.s.  
Excipients : q.s.

### DIRECTION FOR USE:

- For Activation: Pour entire content of activator vial in 5Ltr. Glutaraldehyde solution. Stir well till activated solution colour changes to light blue.
- For Disinfection immerse cleaned medical instruments / equipments completely in activated Glutaraldehyde solution.

### SOAK TIME:

- For high level disinfection (HLD): 20 minutes at 20°C.
- For Sporicidal Activity: 8 – 10 hours at 25°C.
- Wash the instruments thoroughly with sterile water before application.

### AREAS OF APPLICATION:

Glutaraldehyde Solution has excellent material compatibility with heat-resistant & heat sensitive instruments. It is suitable for all common manual and semi-automatic circulation procedures as well as for the fully automatic (Cold Disinfection) procedure. It is also used in ultrasonic baths. It is also suitable for reprocessing of instruments made of rubber, plastic, ABS, metal, porcelain and glass.

**KEEP OUT OF THE REACH OF CHILDREN**

**FOR EXTERNAL USE ONLY**

Storage condition:  
Store in a cool and dry  
place 15°C - 30°C (59°F – 86°F)

### PRECAUTIONS:

Contains Glutaraldehyde.  
Harmful by inhalation and if swallowed. Irritating to respiratory system and skin.  
Risk of serious damage to eyes.  
May cause sensitization by inhalation and skin contact.  
In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Wear suitable PPEs in case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Use only in ventilated areas (refer to the material safety data sheet (MSDS) for additional information). Avoid release to the environment.

### FEATURES:

- High level disinfectant (HLD) for use in reprocessing heat sensitive medical devices.
- Virucidal, Bactericidal, Fungicidal, Tuberculocidal and sporicidal.
- High degree of material compatibility for Medical Devices and Endoscopes.
- Compatible with different water qualities.
- Effective against bacterial spores.
- Can be used for cleaning and disinfection of heat resistant and heat sensitive instruments (Flexible Endoscopes) in the manual treatment, circulation procedures as well as in the Ultrasonic path.



**KERALA GOVERNMENT SUPPLY  
NOT FOR SALE**

Activation Date

Deactivation Date of Activated Solution

The activated solution remains stable  
for 14 days or discard solution if  
colour changes or soiled.



Manufactured in India by  
**Kerala State Drugs and Pharmaceuticals Ltd.**  
Kalavoor P.O., Alappuzha, Kerala-688 522.

Mfg. Lic. No. 9/25/74

Batch No.

Mfg. Date:

Exp. Date:



KERALA STATE DRUGS AND PHARMACEUTICALS LTD.

KALAVOOR PO, ALAPPUZHA, KERALA-688522

Name of the Material: **SPECIFICATION FOR STICKER LABEL**

**GLUTARALDEHYDE SOLUTION I.P. 2% - KMSCL SUPPLY**

SOP No:	KSDP/SOP/SPEC02	Spec. No.:	NB/SPEC/LIQ/SL/018
		Revision No:	01
Effective date:	09/01/2025	Item code	41400071
		Change control No.	25/NBL/003

- Dimension: 170 x100 mm (L X B).
- Base colour: Grey Dawn
- Printing colour: 3 colour
- Sticker label paper: Chromo paper (80GSM).
- Sticker labels should be in sheet form instead of roll(manual labelling).
- Not more than 10 sticker labels in a single sheet
- Sticker label should be perfectly removable from the base paper (half die punch type).
- Sample sheet should be provided for approval.

*Note: Samples to be provided for approval.*

Prepared by	Checked by	Reviewed by	Approved by	Authorised by	<b>MASTER COPY</b> <b>UNCONTROLLED COPY</b>
Officer QA	Production Officer	Production Head	QC Head	QA Head	
Date: 09/01/25	Date: 09/01/2025	Date: 09/01/2025	Date: 09/01/2025	Date: 09/01/2025	