



Kerala State Drugs and Pharmaceuticals Ltd
(A Govt. of Kerala Enterprise)
KALAVOOR, ALAPPUZHA

Phone :0477-2258184 (Extn:207) ; Email :purchase@ksdp.in

Web :http://www.ksdp.co.in

KSDP/PS/EQ/2025-26/400332/T-184

5-2-2026

Notice Inviting Email Quotation for the Supply of Sticker labels and Carton label

Quotations are invited for the supply of undermentioned goods/ services as per the attached specifications on FOR destination basis at our factory site at Kalavoor,Alapuzha, Kerala state

| SI No. | ITEM CODE & DESCRIPTION | UNIT | QUANTITY | EMD | TENDER FEES |
|--------|--|------|-----------|------------------|------------------|
| 1 | 41400071 STICKER LABEL GLUTERALDEHYDE SOLUTION IP 5L CAN IN SHEETS (Blue) | NOS | 5,700 NOS | ♦ Not Applicable | ♦ Not Applicable |
| 2 | 41300082 CARTON LABEL KMSCL GLUTERALDEHYDE (Blue) | NOS | 1,200 NOS | | |
| 3 | 41400070 Sticker Label-Activator for Gluteraldehyde Solution IP 2% -KMSCL(IN SHEETS) (Blue) | NOS | 4,500 NOS | | |

Quotations should be submitted as per the Proforma given below on your letter head.

| SI No. | Name of item | Make | Rate per unit (including freight, if any) | GST % | Offer validity | Remarks, if any |
|--------|--------------|------|---|-------|----------------|-----------------|
|--------|--------------|------|---|-------|----------------|-----------------|

Due Date/Time : 9-Feb-26 / 1 PM

Opening Date/Time: 9-Feb-26 / 1.30 PM

Offer validity*: Minimum 7 days offer validity from the date of closure of bid submission.

Quotation received with offer validity less than 7 days from the date of closure of bid submission will be entirely rejected.

Please send your lowest offers of the item to our e-mail ksdptender@gmail.com before 1 PM, 09/02 /2026.

The Quotation should be submitted through a password protected excel sheet. Please share your Password to our email ksdptender@gmail .com@1.30 PM, 09.02.2026.

TERMS & CONDITIONS

- 1.Payment terms :- 30 days after the receipt of the material along with documents, subject to QC approval.
- 2.Mode of payment:- E-Payment.
- 3.Delivery Period:- Supply should be effected within 15 days on award of PO
- 4.Special instructions:- Quotation number Should be mentioned in the subject line.

HOD - Purchase



KERALA STATE DRUGS AND PHARMACEUTICALS LTD.

KALAVOOR PO, ALAPPUZHA, KERALA-688522

Name of the Material: **SPECIFICATION FOR STICKER LABEL**

GLUTARALDEHYDE SOLUTION I.P. 2% - KMSCL SUPPLY

| | | | |
|-----------------|-----------------|--------------------|--------------------|
| SOP No: | KSDP/SOP/SPEC02 | Spec. No.: | NB/SPEC/LIQ/SL/018 |
| | | Revision No: | 01 |
| Effective date: | 09/01/2025 | Item code | 41400071 |
| | | Change control No. | 25/NBL/003 |

- Dimension: 170 x100 mm (L X B).
- Base colour: Grey Dawn
- Printing colour: 3 colour
- Sticker label paper: Chromo paper (80GSM).
- Sticker labels should be in sheet form instead of roll(manual labelling).
- Not more than 10 sticker labels in a single sheet
- Sticker label should be perfectly removable from the base paper (half die punch type).
- Sample sheet should be provided for approval.

Note: Samples to be provided for approval.

| Prepared by | Checked by | Reviewed by | Approved by | Authorised by | MASTER COPY UNCONTROLLED COPY |
|----------------|--------------------|------------------|------------------|------------------|--|
| Officer QA | Production Officer | Production Head | QC Head | QA Head | |
| | | | | | |
| Date: 09/01/25 | Date: 09/01/2025 | Date: 09/01/2025 | Date: 09/01/2025 | Date: 09/01/2025 | |

170 mm

100 mm

5 Litre

GLUTARALDEHYDE SOLUTION I.P. 2%

Composition: Glutaraldehyde 2%w/v
Purified Water I.P. : q.s.
Excipients : q.s.

DIRECTION FOR USE:

- For Activation: Pour entire content of activator vial in 5Ltr. Glutaraldehyde solution. Stir well till activated solution colour changes to light blue.
- For Disinfection immerse cleaned medical instruments / equipments completely in activated Glutaraldehyde solution.

SOAK TIME:

- For high level disinfection (HLD): 20 minutes at 20°C.
- For Sporicidal Activity: 8 – 10 hours at 25°C.
- Wash the instruments thoroughly with sterile water before application.

AREAS OF APPLICATION:

Glutaraldehyde Solution has excellent material compatibility with heat-resistant & heat sensitive instruments. It is suitable for all common manual and semi-automatic circulation procedures as well as for the fully automatic (Cold Disinfection) procedure. It is also used in ultrasonic baths. It is also suitable for reprocessing of instruments made of rubber, plastic, ABS, metal, porcelain and glass.

KEEP OUT OF THE REACH OF CHILDREN

FOR EXTERNAL USE ONLY

Storage condition:
Store in a cool and dry
place 15°C - 30°C (59°F – 86°F)

PRECAUTIONS:

Contains Glutaraldehyde.
Harmful by inhalation and if swallowed. Irritating to respiratory system and skin.
Risk of serious damage to eyes.
May cause sensitization by inhalation and skin contact.
In case of contact with eyes, rinse immediately with plenty of water and seek medical advice.

Wear suitable PPEs in case of accident or if you feel unwell, seek medical advice immediately (show the label where possible).

Use only in ventilated areas (refer to the material safety data sheet (MSDS) for additional information). Avoid release to the environment.

FEATURES:

- High level disinfectant (HLD) for use in reprocessing heat sensitive medical devices.
- Virucidal, Bactericidal, Fungicidal, Tuberculocidal and sporicidal.
- High degree of material compatibility for Medical Devices and Endoscopes.
- Compatible with different water qualities.
- Effective against bacterial spores.
- Can be used for cleaning and disinfection of heat resistant and heat sensitive instruments (Flexible Endoscopes) in the manual treatment, circulation procedures as well as in the Ultrasonic path.



**KERALA GOVERNMENT SUPPLY
NOT FOR SALE**

Activation Date

Deactivation Date of Activated Solution

The activated solution remains stable
for 14 days or discard solution if
colour changes or soiled.



Manufactured in India by
Kerala State Drugs and Pharmaceuticals Ltd.
Kalavoor P.O., Alappuzha, Kerala-688 522.

Mfg. Lic. No. 9/25/74

Batch No.

Mfg. Date:

Exp. Date:



KERALA STATE DRUGS AND PHARMACEUTICALS LTD.

KALAVOOR PO, ALAPPUZHA, KERALA-688522

PACKING SPECIFICATION

LABEL SPECIFICATION(MASTER CARTON LABEL) – KMSCL SUPPLY

| | | | |
|-----------------|-----------------|--------------|-----------------------|
| SOP No: | KSDP/SOP/SPEC02 | Spec. No: | KSDP/SPEC/KMSCL/CL/68 |
| Effective date: | 25/05/2023 | Revision No: | 00 |

| Sl. No. | NAME OF THE PRODUCT | MASTER CARTON LABEL | | |
|---------|---------------------------------------|---------------------|--------|-----------|
| | | DIMENSION | COLOUR | PAPER GSM |
| 1 | GLUTARALDEHYDE SOLUTION I.P.2% 5L CAN | 15 x 10 cm | 2 | 60 |

- Maplitho paper
- Labels are in single piece, not in roll or sheet form.

Samples should be provided for approval.

| Prepared by | Checked by | Reviewed by | Approved by | Authorised by | MASTER COPY UNCONTROLLED COPY |
|----------------|--------------------|-----------------|----------------|----------------|----------------------------------|
| Officer QA | Production Officer | Production Head | QC Head | QA Head | |
| | | | | | |
| Date: 20/05/23 | Date: 23/05/23 | Date: 25/05/23 | Date: 25/05/23 | Date: 25/05/23 | |

GLUTARALDEHYDE SOLUTION I.P. 2%

ഗ്ലൂട്ടറാൾഡിഹൈഡ് സൊല്യൂഷൻ ഐ.പി. 2%

5Litre

Mfg. Lic. No.: 9/25/74

Batch No. :

Mfg. Date :

Exp. Date :

Qty. Packed:

Net Weight :kg.



കേരള സർക്കാർ വിതരണം - വില്പനയ്ക്കുള്ളതല്ല



Manufactured in India by
Kerala State Drugs and
Pharmaceuticals Ltd.

Kalavoor P.O., Alappuzha
Kerala - 688 522

KERALA GOVT. SUPPLY
NOT FOR SALE



KERALA STATE DRUGS AND PHARMACEUTICALS LTD.

KALAVOOR PO, ALAPPUZHA, KERALA-688522

Name of the Material: **SPECIFICATION FOR STICKER LABEL**

ACTIVATOR FOR GLUTARALDEHYDE SOLUTION I.P. 2% - KMSCL SUPPLY

| | | | |
|-----------------|-----------------|--------------|----------------------|
| SOP No.: | KSDP/SOP/SPEC02 | Spec. No.: | NB/SPEC/LIQ/SL/018/A |
| Effective date: | 18/07/23 | Revision No: | 00 |

- Dimension: 105 x30 mm (L X B).
- Base colour: White
- Printing colour: Black.
- Sticker label paper: Chromo paper (80GSM).
- Sticker labels should be in sheet form instead of roll(manual labelling).
- Sticker label should be perfectly removable from the base paper (half die punch type).
- Sample sheet should be provided for approval.

Note: Samples to be provided for approval.

| Prepared by | Checked by | Reviewed by | Approved by | Authorised by | MASTER COPY UNCONTROLLED COPY |
|----------------|--------------------|-----------------|----------------|----------------|----------------------------------|
| Officer QA | Production Officer | Production Head | QC Head | QA Head | |
| | | | | | |
| Date: 10/07/23 | Date: 11/07/23 | Date: 12/07/23 | Date: 12/07/23 | Date: 13/07/23 | |

Activator for Glutaraldehyde Solution I.P. 2%

Direction: Put the content of the bottle in the 5 litre Glutaraldehyde solution I.P. 2 % and shake the solution well to dissolve completely.